

	National Institute for Occupational Safety and Health National Personal Protective Technology Laboratory Evaluation and Testing Branch
---	--

TEST REPORT

Task Number: 27899

Manufacturer: Rizhao Sanqi Medical & Health Articles Co., Ltd.

Prepared by: Bryan Tennant

Date: September 16, 2025

Tests Conducted by: Matt Duling and Bryan Tennant

Respirator Tested: SQ100D

Tests Completed

<u>Test Description</u>	<u>STP Number</u>
A. Exhalation Resistance Test	TEB-APR-STP-0003
B. Inhalation Resistance Test	TEB-APR-STP-0007
C. Sodium Chloride (NACL) N-95 Test	TEB-APR-STP-0059

Overall Results

The respirator system tested met the requirements of all the above procedures.